

# STAINLESS STEEL FURNITURE & FUME HOODS



## Quality Assured

Mott's stainless steel product offering is designed and manufactured to deliver practical solutions to specialized lab environments.

Our entire painted steel product line is available in type 304 or 316 stainless steel and is manufactured with a No.4 finish. Type 316 is made with Molybdenum to provide an extra level of corrosion resistance.

Applications include pharmaceutical manufacturing, research labs, clean rooms, food processing, health care, operating rooms, animal facilities, autopsy areas, and radioactive processes.

## Moisture & Corrosion Resistance

Properly selected and maintained stainless steel provides excellent moisture related corrosion resistance.

## Impact Resistance

Stainless steel's high yield strength, even at elevated temperatures, provides high impact strength and therefore ensures durability and extended product life.

## Sterile & Easy to Sanitize

Stainless steel work surfaces with smooth radius corners are ideal for environments requiring high levels of cleanliness.



## Attractive & Aesthetic Appearance

Stainless steel will compliment your lab's design or color scheme and provide a cohesive looking laboratory. It reflects light, offering dark labs an easy-to-work-in setting. The easily maintained surface of stainless steel provides an attractive and timeless appearance.

## Friendly to the Environment

Stainless steel is 100% recyclable and has value even after a very long life in the laboratory. Approximately 75% - 85% of new stainless steel comes from recycled stainless steel.



# Resilient

If your lab demands a hygienic and sterile space, then stainless steel is a top choice. Stainless steel is durable, long-lasting and easy to sanitize due to its non-porous qualities. It is also resistant to corrosion and heat, making it an excellent choice for laboratory countertops and sinks.



## Casework & Tables

Built to withstand challenging laboratory requirements, Mott offers a wide variety of fixed and flexible solutions in many configurations with a wide variety of options. Our casework and furniture will enhance any laboratory environment and meet your specific aesthetic needs.



## Fume Hoods

All Mott's fume hoods are readily available in stainless steel construction. Features can include coved corner liners with integral work surfaces where high levels of cleanliness are required. Hoods can be custom designed and manufactured to meet specific application requirements.



## Worksurfaces, Scullery & Scrub Sinks

Built to your exact requirements, worksurfaces and sinks can be configured with features such as single or multiple compartments, drain boards, side, and back splashes, sensor faucets, foot pedal valves and more...



## Customized Solutions

Mott offers a variety of specialty products to fit the needs of most labs, such as: operating room casework, nurse's stations, biological safety casework, and more. Mott is always ready to design and manufacture custom engineered products.