

OPTIMA™ MOBILE WORKSTATION

Flexible Solutions

Mobile Optima™ self contained workstations are designed to enhance lab productivity. Push a button and the work surface adjusts to meet the height demands of individual users and equipment. Add plumbing, data, shelving, and accessories to meet your specific requirements. You can easily reconfigure lab layouts as needs change; simply move it where you need it when you need to.

Configurations - Assembly is 84" high with 15" deep uprights. Available in single sided and double sided assemblies, in various widths.

Illuminate The Workspace - Standard integrated overhead task lights delivers bright light to the workspace.

Integrated Services - Light switch, duplex receptable, and worksurface height button located on the inside face of the upright.

Push-Button Adjustable Worksurface - Raise or lower the worksurface from 30" to 36" with a touch of a button with the electric height adjustable table frame.

Additional Electrical & Data - Two duplex receptacles are provided on the lower service chase of the single-sided assembly. Double-sided assembly receives two duplex receptacles on each side. Data ports are an available option on the lower service chase.

Service Accessibility - Hinged panel on the service chase is easily accessed for convenient service maintenance (located on other side of assembly).



Single-sided assembly shown

Overhead Service Distribution - Workstation provides quick and easy connection to Mott's overhead service solutions.

Upper Storage - Optional fixed shelves mount between the uprights to optimize vertical space and allow for maximum utilization of valuable work surface space.

Accessories - Optima™ Workstations can be configured with many accessories such as keyboard trays, monitor arms, and shelf task lights to suit your needs.

Plumbing Services - Up to three factory pre-plumbed service fixtures available. The fixtures are easily accessible while keeping the work surface open.



Removable End Panels - Designed for rugged durability, these panels provide an inspection point for connections, shut-offs and service lines. Access for modifications is quick and easy.

Effortless Mobility - Swivel casters have an adjustment knob that lowers the pad and raises the wheel off the ground to prevent rolling. Perfect for applications where work demands stability and a level work surface.



mott
GROUP OF COMPANIES

Steel Division & HQ - 452 Hardy Rd, Brantford ON, Canada N3T 5L8 • T. 519.752.7825
Panel Products Division - 562 Industrial Park Rd, Maxwelton WV, USA 24957 • T. 304.497.2115

www.mott.ca

QUALITY BY DESIGN