

LABORATORY TABLES

Versatile. Multi-Functional.



Looking For The Perfect Laboratory Table?

Whether it's for a teaching lab or a cutting-edge research facility, Mott has a wide range of adaptive solutions for creating flexible lab spaces to meet present and future demands. Mott has a solution for you.



Versatile

Mott's collection of tables are designed to support your lab and the way you work:

- basic lab tables
- adjustable height tables
- height adjustable hydraulic tables
- all-welded steel tables
- heavy-duty tables
- mobile tables
- research and educational tables

Combine numerous options, in a range of materials, and make a table that fits your needs:

- pin or bolt height adjustable legs
- various types of casters
- easy glide nylon levelers
- electrical
- vibration dampening or vibration isolation
- keyboard trays and numerous storage options
- hanging rails for suspended casework
- fixed shelves



Wood Sturdy & Warm

We manufacture to the highest standards using architect selected materials on every project and can meet your appearance, durability, and environmental needs.



Stainless Steel Sanitized

If your lab demands a hygienic and sterile space, then stainless steel is a top choice. Stainless steel is durable, long-lasting, and easy to sanitize due to its non-porous qualities. It is also resistant to some corrosives and heat.



Collaboration

The best lab spaces support collaborative and individual work environments. The Mott collection of tables creates a range of adaptable and high performing space for users and tasks across the entire lab.



Research & Education

We can design custom tables to meet the special needs of your research or teaching lab. Tables are available with integrated fixtures and accessories, as well as cabinets to store your lab equipment in one space-saving location.



Mobility

Based on the mobility needs of the user, most tables can be provided with various casters, providing effortless mobility around the lab. Optional heavy-duty leveling casters have an adjustment knob for times when you need both mobility and absolute stability.



Fully Welded

When extra strength is required to support heavy loads and lessen vibration, choose from a selection of all-welded, height adjustable table frames in powder coated steel and stainless steel; add casters for more mobility.



Height Adjustable

Table frames can be quickly adjusted in 1” increments to easily adapt to new equipment and/or users. Tables are available in powder coated steel, stainless steel, and wood.

Instantly Adjustable

Raise or lower the work surface at the touch of a button with an electric height adjustable option. Or choose the crank-operated table frame option; legs are supplied with leveling feet. Suspended cabinets or shelving are also available options.



Sharing Resources & Maximizing Utilization

Relocatable tables and cores carrying multiple services, all within a sturdy UL listed system.

Affordable

Tables are a cost-effective furniture system that lab personnel can reconfigure quickly and efficiently.



Need Extra Storage?

Most tables can accommodate our full range of suspended cabinets, drawers, and shelves. Other lower storage options include mobile cabinets and floor mounted cabinets on glides.



Flexible

We offer tables with endless flexibility to meet your changing needs. Reconfigure quickly and easily.



Craftsmanship

Mott is an AWI QCP approved manufacturer and our cabinets are built and finished to SEFA standards.

Sustainability

Offering FSC certified and sustainable materials and no added urea formaldehyde (NAUF) products. We are California Air Resources Board (CARB) compliant.



Service Distribution

Overhead service carriers and ceiling panels are designed to complement moveable furniture systems by providing a convenient point of access for mechanical, electrical, data, and plumbing connections at the ceiling. Available in a variety of sizes, types and styles, there is a carrier to meet your specific lab requirements.

Products **designed** for your lab.
Explore our **solutions**.

[**www.mott.ca**](http://www.mott.ca)