

GLASSGUARD™

TABLE TOP FUME HOOD

The GlassGuard™ table top fume hood is a compact, easy to install solution for less demanding laboratory activities which require ventilation. It may be used for a variety of smaller scale activities including weighing and dispensing, package opening and various instructional activities. It is connected to the laboratory exhaust system. Exhaust air required is considerably less than a full size chemical fume hood.

Exhaust Collar - 6" round plastic exhaust duct for connection to laboratory exhaust system or remote fan (not included).

Duct Connection - Exhaust duct connection may be changed on site to be exhausted up or down.

Cord Pass-Through - Two utility ports allows equipment to be plugged into receptacles on adjacent benchtops.

Interior Workspace - Usable interior workspace depth is 23". Work top not included.

Airfoil - Compact raised airfoil across front of countertop directs air rearward.

Configurations - Assembly is 32" high by 30" deep. Available in sizes 24", 30", 36", and 48" wide.



Illuminate The Workspace - Optional LED lighting powered by external low voltage power supply, delivers bright light to the workspace.

Laminated Safety Glass - Sides, top and sash provide a well-lit workspace with excellent visibility and impact resistance.

Sash Opening - When the sash is down for operation, the opening is 8.75".

User Friendly - Ergonomic pivoting sash is sloped 20° to allow user to overlook the working area.



Pivoting Sash - The sash pivots upward to allow tall apparatus to access work chamber providing a vertical opening height of 19".

Coating Finish and Color - A thermosetting laboratory grade powder coating delivers exceptional chemical and scratch resistance, excellent hardness, adhesion, and a quality luster. The long lasting finish is available with 15 standard color options, including multiple color schemes, and color matching capabilities.

Agency Approvals - Tested to ASHRAE 110-2016 and meets industry standards for fume hood containment.



Base Configurations - Place on any suitable sized flat benchtop; we offer a full line of cabinets and tables to compliment the table top hood.