

# FURNITURE & FUME HOODS FOR EVERY PURPOSE

Designed For **Labs**  
Designed For **people**





**At Mott, we believe our products  
make a difference in your lab.**

That's why we make furniture and fume hood solutions that work for you.

## Who We Are

Mott places great importance on the design of our products, the health and well-being of our customers, our employees, our community and the environment. Delivering innovative ways to improve the performance of our customer's laboratories has become our mission.

## 60 Years Of Innovative Design And Manufacturing

Mott's laboratory furniture and fume hoods can be found in tens of thousands of laboratories in educational, pharmaceutical, biotech, health care, chemical, and industrial laboratories across North America and the world.



## Built To Last

### World-Class Manufacturing

Excellence and innovation is built into every design. Our ISO system ensures there is quality in everything we do. Mott is committed to technology. Each of our manufacturing facilities are outfitted with the highest state-of-the-art tooling technology, allowing us to manufacture quality products on time and affordably. Our products have been tested to meet SEFA Recommended Practices to ensure many years of intensive use in a demanding environment and we're continually improving our designs and processes.



### Designing For Sustainability

Mott is committed to protection of the environment, conservation of natural resources, and energy. We're focusing on corporate environmental responsibility while enhancing product quality with modern, efficient, and environmentally responsible manufacturing.

Our objective is to future-proof the laboratory for changes that are anticipated or unanticipated. Our products and systems are designed to be reconfigured and relocated, offering flexibility to labs. Our products are durable, easy to maintain, and we supply products compatible with those purchased years ago. This all supports operational functionality for an extended life cycle.







## Quality **Craftsmanship**

Mott is quality certified by the Architectural Woodworking Institutes (AWI) highest grade of certification. FSC certified and no added urea-formaldehyde (NAUF) are available.



# Lab Casework Solutions

Mott's casework is designed and manufactured to provide durable solutions to meet the unique requirements of your laboratory.

Casework can be provided in a variety of materials:

- Powder Coated Steel
- Powder Coated Steel with Wood Veneer, Plastic Laminate, or Phenolic Resin Fronts
- Stainless Steel
- Wood Veneer
- Plastic Laminate
- Phenolic Resin



## Tailored Solutions

Available in inset or overlay design with numerous configurations and countless options, Mott's casework can be configured for fixed, mobile, or suspended installations. Your options for steel and wood include a wide variety of standard colors and stains, as well as the ability to custom match. A wide spectrum of choices are available for plastic laminate and phenolic resin.



## Enhanced **Protection**

Our steel casework can be constructed from satin galvanized steel for an extra level of protection against rust and corrosion in damp or humid environments.





## Workspace **Adaptability**

Integrate plumbing, electrical and data for quick and easy access to services. Systems can be configured with a broad range of accessories, suspended cabinets, task lighting, pegboards, marker boards and tack boards.



## **Flexible & Mobile Systems**

We offer a wide range of adaptive solutions for creating flexible lab spaces to meet present and future demands. Systems that enhance lab productivity and space efficiency; Mott has a solution for you.

Choose the system that best suits your needs:

- Moveable Benches
- Mobile Carts & Workstations
- Tables
- Core Systems
- Service Distribution Carriers



## **Easily Reconfigured**

Flexible and moveable systems allow labs to be easily rearranged and increases mobility within the lab. Our systems work together, providing a great blend of possibilities to your laboratory layout.



## **Ergonomics**

Lab users deserve productive, comfortable, and safe workspaces. Mott's flexible and adaptable furniture systems are designed from the ground up to promote productivity and enhance worker ergonomics. Lateral and vertical adjustability allows workers the freedom to create work areas that meets the needs of their processes while maximizing their personal comfort and safety.





## Fume Hood Specialists

We make thousands of fume hoods every year with a very broad selection of cataloged hoods and options:

- General Chemistry
- High Performance
- Specialty Applications
- Teaching
- Tabletop
- Automatic Sash Operators
- Integrated Variable Air Volume Systems
- Canopy Hoods
- Ventilated Workstations

No manufacturer offers more **customization** with less hassle.



## Safety

Our strong commitment to fume hood safety, coupled with unparalleled technical expertise, serves as the foundation for our long-standing customer relationships. Our hoods meet these stringent performance requirements: UL Classified to UL1805, CSA Certified to CSA C22.2 No. 61010, UL Standard 61010-1, and tested in accordance with ASHRAE 110.

We also offer standard, third party, and specialty investigative testing at our on-site fume hood test facility.

Mott builds fume hoods that are performance ready for your tough requirements.

## We Focus On Your Needs

Our experts have a vast knowledge of designing safe, durable and energy efficient hoods for almost every application. Our attention to detail has allowed Mott's fume hoods to earn a reputation as an industry leader in reliable containment performance. Our specialists will evaluate your needs and recommend the right fume hoods for your lab.







**Yes we can!**  
**Customization, it's our specialty.**

What's your vision for tomorrow? No matter the task, we can design and build it. We have a proven track record in the market as the company who can turn an idea into a full solution, and quickly.



## Work Surfaces, Sinks & Pegboards

Mott's extensive stainless steel experience can deliver practical solutions to your specialized laboratory environments.

Custom designed and manufactured to your exact requirements:

- Stainless Steel Work Surfaces
- Phenolic Resin Work Surfaces
- Stainless Steel Scullery Sinks
- Pegboards



## Resilient

If your lab demands a hygienic and sterile space, then stainless steel is a top choice. Stainless steel is durable, long-lasting and easy to sanitize due to its non-porous qualities. It is also resistant to corrosion and heat, making it an excellent choice for laboratory countertops and sinks.



**When you need experience, our dealers have it covered.**

Mott has formed a strong network of dealers who are knowledgeable, reliable, and flexible, with a clear focus on customer needs. They have decades of experience in drafting, designing, layout, budgeting and estimating, procurement, project planning, installations, and much more.

## Trust Us

It's not just words; it's one of our company core values. At Mott we take great pride in our high standards of customer care, after all, customers are our greatest assets. So we feel it is important to ensure customer satisfaction is our number one priority.





Products **designed** for your lab.  
Explore our **solutions**.

[\*\*www.mott.ca\*\*](http://www.mott.ca)