

# GLASSGUARD™ VENTILATED WORKSTATION

The GlassGuard™ table top ventilated workstation is a compact, easy to install solution for less demanding laboratory activities which require ventilation. It may be used for a variety of smaller scale activities including weighing and dispensing, package opening and various instructional activities. It is connected to the laboratory exhaust system. Exhaust air required is considerably less than a full size chemical fume hood.

**Exhaust Collar** - 6" round plastic exhaust duct for connection to laboratory exhaust system or remote fan (not included).

**Duct Connection** - Exhaust duct connection may be changed on site to be exhausted up or down.

**Cord Pass-Through** - Two utility ports allows equipment to be plugged into receptacles on adjacent benchtops.

**Interior Workspace** - Usable interior workspace depth is 23". Work top not included.

**Airfoil** - Compact raised airfoil across front of countertop directs air rearward.

**Configurations** - Assembly is 32" high by 30" deep. Available in sizes 24", 30", 36", and 48" wide.



**Illuminate The Workspace** - Optional LED lighting powered by external low voltage power supply, delivers bright light to the workspace.

**Laminated Safety Glass** - Sides, top and sash provide a well-lit workspace with excellent visibility and impact resistance.

**Sash Opening** - When the sash is down for operation, the opening is 8.75".

**User Friendly** - Ergonomic pivoting sash is sloped 20° to allow user to overlook the working area.



**Pivoting Sash** - The sash pivots upward to allow tall apparatus to access work chamber providing a vertical opening height of 19".

**Coating Finish and Color** - A thermosetting laboratory grade powder coating delivers exceptional chemical and scratch resistance, excellent hardness, adhesion, and a quality luster. The long lasting finish is available with 15 standard color options, including multiple color schemes, and color matching capabilities.

**Performance** - Factory tested to ASHRAE 110-2016 and exceeds industry standards for containment.



**Base Configurations** - Place on any suitable sized flat benchtop; we offer a full line of cabinets and tables to compliment the table top hood.